



Daves mobility 336-534-7722

Q6 edge[®]
HD



Get the Edge HD with iLevel!

With standard features like 4-pole motors, Mid-Wheel 6[®] Design, and ATX Suspension, and optional iLevel[®] power adjustable seat height, the Q6 Edge[®] HD Power Chair is engineered to meet the performance needs of the most active

user. The Q6 Edge HD Power Wheelchair accepts a complete range of seating and electronics options to deliver exceptional rehab capability at an excellent value.

Q6 Edge HD Features

- Mid-Wheel 6® allows six wheels on the ground for maximum stability
- Compatible with TRU-Balance® 3 HD Power Positioning Systems
- Available with Q-Logic 3 Advanced Drive Control System, which offers Bluetooth® access to computer, tablet and phone functions
- 450 lbs. weight capacity
- Low impact OMNI-Casters (nylon, spherical-shaped casters) on front and rear prevent wheel hang-ups
- ATX Suspension (Active-Trac® with extra stability) incorporates front OMNI-Casters and semi-independent rear caster beam for enhanced performance over more varied terrain
- 400 lbs. weight capacity iLevel® option enables up to 10” of seat elevation while at walking speed (up to 3.5 mph).
- Optional iLevel's patented Extra Stability Technology® enhances safety while performing transfers, reach, and numerous other activities of daily living
- Optional LED fender lights (standard with iLevel)
- Available with WC19 crash-test approved occupied transit kit option on Synergy® seating



NE Controller

- Non-expandable rehab electronics deliver simplicity without limiting features

- Features an easy-to-read data display, soft-touch buttons and charger port
- For use with power chairs without power seating options

• Q6 Edge HD Specifications

Model	Q6 Edge HD
Drive Wheels	14" run flat Optional: 14" Pneumatic, 14" Knobby
Caster Wheels	Front: 6" Solid Rear: 6" Solid
Anti-Tip Wheels	N/A
Maximum Speed⁴	Up to 4.7 mph
Ground Clearance	3" (frame center)
Turning Radius⁵	19.75"
Overall Length	35.6" without foot riggings
Base Width	26.25"
Seating Sizes	TRU-Balance® 3 HD Synergy®/Static: W: 20-32" D: 20-28" Power Tilt: W: 20-32" D: 20-25" Lift & Tilt: W: 20-24" D: 20-25"

Seat-to-Floor Heights¹⁶	TRU-Balance® 3 HD Synergy®/Static ⁸ : 17.5-19.5" Power Tilt ⁸ : 18.5-20.5" Lift & Tilt ⁸ : 17.75-19.75"
Manual Recline¹⁰	No
Battery Size⁶	Group 24
Battery Weight¹⁴	52 lbs.
Available Electronics¹²	75A Q-Logic 2 NE 75A Q-Logic 2 NE+ 120A Q-Logic 3 EX
Battery Charger	8A Off-board
Motor Packages	4-Pole 4.7 mph ¹
Weight Capacity⁷	450 lbs. ^{15, 16}
Base Weight	143.5 lbs.
Medicare Code Coverage	Group 3: Heavy Duty (K0850/K0851) Single Power HD (K0858) Multiple Power HD (K0862)

Notes:

¹ Standard

⁴ Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.

⁵ Without front riggings

⁶ Battery capacity may change when adding power positioning systems.

⁷ Weight capacity is dictated by the selection of seating.

⁸ All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Depths may increase seat-to-floor height. Contact Customer Service for details. Limited measurements are in the lowest position; does not include elevating system unless specified.

¹⁰ Maximum dimensions of the Manual Recline vary by power base and seating system. Call for details.

¹² All Quantum Rehab® units can be ordered with specialty controls

¹⁴ Battery weight may vary +/- 3 lbs. based on manufacturer.

¹⁵ Weight capacity is reduced to 400 lbs. with the iLevel® power adjustable seat height.

¹⁶ A minimum seating size of 20" x 20" is required with iLevel®. iLevel maximum seating size is 24" x 25".

¹⁶ Seat-to-floor heights are measured with the unit unloaded.